

NEW




Stuart®

---

## Mini-Stirrer

Stuart® would like to introduce its new mini stirrer. The SM5/ stirrer compliments its comprehensive range of hotplate stirrers, offering a fun and practical powerful stirrer. The mini stirrer is completely new and is available in a range of designs. Capable of stirring up to 1 litre, the stirring speed is adjustable up to 2,000rpm.



## Stuart Mini Stirrer, SM5/



- Lightweight stirrer
- Capable of stirring up to 1L
- Universal power adaptor
- Choice of designs
- With BioCote antimicrobial protection

The SM5/ mini stirrer from Stuart is a powerful magnetic stirrer, capable of stirring up to 1 litre, the stirring speed is adjustable up to 2,000rpm. The mini stirrer is available in a choice of designs, please ensure you use the correct order code to receive your chosen model. The unit is made from robust polypropylene with a chemically resistant polycarbonate top.

## Technical Specification

Capacity, Litres	1
Speed, rpm	350 to 2,000
Weight, g	500
Dimensions, w x d x h, (mm)	143 x 143 x 66
Electrical supply	120-230V, 50-60Hz
IP Rating	IPX1

## Ordering Information

<b>Model</b>	<b>Description</b>
SM5/STUART	Stuart mini stirrer, Stuart fascia
SM5/BIBBY	Stuart mini stirrer, Bibby fascia
SM5/SWIRL	Stuart mini stirrer, Swirl fascia

## Contact

Bibby Scientific - UK (Group HQ)  
Beacon Road, Stone, Staffordshire,  
ST15 0SA, United Kingdom  
Tel: +44 (0)1785 812121  
Fax: +44 (0)1785 813748  
e-mail: [sales@bibby-scientific.com](mailto:sales@bibby-scientific.com)

[sales@bibby-scientific.com](mailto:sales@bibby-scientific.com)



**Find out more!!!**

Please scan the QR/Mobile  
Tag with your smartphone  
for more information